A new standard. A new choice.

Philips Respironics Nuance Pro gel pillows mask

Mask fitting guide



Before putting on the mask, place the pillow cushion tips under your nostrils. Choose the cushion size that sits comfortably under your nose without any gaps.

Important: Do not push the pillows gel base into your nostrils.

Pulling the headgear over your head, gently hold the pillows cushion under your nose. The headgear should lie flat against the back of your head.





Position the pillow cushion tips into your nostrils so that they sit comfortably under your nose without any gaps.

fit loosely and feel comfortable.





Adjust the side straps evenly by pulling upward; the gel pads should sit under your cheek bones. The mask should feel secure and comfortable against your face. Do not overtighten.



Final adjustments for comfort and seal should be made while lying down with your device pressure applied.

Note: if excessive leaks persist, try an alternate cushion size.

Additional details on fitting may be found in the instructions that accompany your mask.



For more information, please visit www.philips.com/respironics

Broudy AE 1/14/14 MCI 4105968 PN 1114507

A new standard. A new choice.

Philips Respironics Nuance gel pillows mask

Mask fitting guide



Before putting on the mask, place the pillow cushion tips under your nostrils. Choose the cushion size that sits comfortably under your nose without any gaps.

Important: Do not push the pillows gel base into your nostrils.

Pulling the headgear over your head, gently hold the pillows cushion under your nose. The headgear should lie flat against the back of your head.



Position the pillow cushion tips into your nostrils so that they sit comfortably under your nose without any gaps.



4

Adjust the top crown strap so that the headgear sits just above your ears and is not too close to your eyes. The headgear should fit loosely and feel comfortable.



Adjust the side straps evenly by pulling downwards; the mask should feel secure and comfortable against your face. Do not overtighten.



6

Final adjustments for comfort and seal should be made while lying down with your device pressure applied.

Note: if excessive leaks persist, try an alternate cushion size.

Additional details on fitting may be found in the instructions that accompany your mask.

For more information, please visit www.philips.com/respironics
Broudy AE 1/14/14 MCI 4105968 PN 1114507

PHILIPS

RESPIRONICS

A smaller fabric frame is also available by request.