

# SoClean. & SoClean 2 Go Compatibilty Chart

Most CPAP machines with a standard hose will work with the SoClean right out of the box. Some CPAP machines may need a specific adapter, particularly units with a heated hose. To see whether you need to have a special adapter to use the SoClean with your CPAP, please refer to our guide below.



N/A = Not Applicable    ✓ = Out of Box    + = Need Adapter    X = Not Compatible

		PAP	WITH HUMIDIFIER	WITH HEATED HOSE
<b>Aeiomed</b>				
	Everest	✓	✓	N/A
<b>Apex Medical</b>				
	ICH	✓	✓	N/A
	XT	✓	✓	N/A
<b>Carefusion</b>				
	PureSom	✓	✓	N/A
<b>Devilbiss</b>				
	ComfortPAP	✓	✓	N/A
<b>Fisher &amp; Paykel</b>				
	Icon	✓	✓	+
	SleepStyle 200 Series	✓	✓	N/A
	SleepStyle 600 Series	+	+	+
<b>Philips Respironics</b>				
	PR System One REMstar 60 Series	✓	✓	+
	PR System One REMstar	✓	✓	N/A
	SleepEasy	✓	✓	N/A
	SleepEasy II	✓	✓	N/A
	M-Series V 2	✓	✓	N/A
	M-Series V 3	✓	✓	N/A
	Remstar	✓	✓	N/A
<b>PMI Probasic</b>				
	Zzz-PAP	✓	✓	N/A
<b>Puritan Bennett</b>				
	GoodKnight 420	✓	✓	N/A
	Sandman	✓	✓	N/A
<b>RESmart</b>				
	600	✓	✓	N/A
	700	✓	✓	N/A
<b>ResMed</b>				
	S8	✓	✓	N/A
	S9	✓	✓	+
	AirSense 10	+	+	+
<b>Transcend</b>				
	Transcend II	✓	(HME) X	N/A