

Quattro™ Air ...  
the lighter way  
to better sleep.



ResMed's new lightweight **Quattro Air** Full Face Mask builds on Quattro's proven performance and reliability. It provides enhanced comfort and ease of use, and is 45% lighter\* than the previous Mirage Quattro mask.

## Benefits for you ...

### **Builds on the strength of Quattro**

Quattro Air provides performance you can rely on, in a new sleeker, lighter mask with less bulk on the face, so it's bound to have strong patient appeal.

### **Less parts for quick and easy set-ups**

Quattro Air has fewer parts than most other Full Face Masks (FFM) – just four. It's quick to demonstrate and easy to train patients to use, supporting patient compliance.

### **Improved patient usability**

The mask features greatly reduced complexity with improved comfort and stability, reducing patient call-backs and freeing up more of your time.

### **Innovative design for easy acceptance**

It's light, streamlined and comfortable, allowing you to assist those of your patients who struggle to stick with therapy.

## ... and your patients

### **Lightweight and less bulk on the face**

Quattro Air is light with a clearer field of vision and less bulk on the face than other FFM's, so patients are less likely to feel claustrophobic or intimidated.

### **Intuitive and uncomplicated**

Quattro Air is simple for patients to use, for quick and easy cleaning and reassembly.

### **Forehead stability**

Its advanced Flex-wing forehead support is comfortable and feels light on the face.

### **Quiet operation**

The circular venting gently and quietly diffuses airflow away for less user and partner disturbance—only 28 dBA.\*\*

\* The Medium size of Quattro Air is 45% lighter than the Medium size of Mirage Quattro.

\*\* Declared dual numbered noise emission values (with uncertainty 3 dBA) in accordance with ISO 4871.e1.



**Quattro™ Air**  
FULL FACE MASK

## Even more comfort ...

Quattro Air is lighter, sleeker, and even more comfortable than the previous Mirage Quattro mask due to a series of design enhancements.

It's also simpler to use, trimming down the number of parts to just four. That makes it easier for patients to understand how to use their Quattro Air and simple for them to put it back together after cleaning.

### SoftEdge™ headgear design

has premium rolled-edge fabric that reduces facial marking and is also comfortable to wear



### Improved comfort

and less chance for skin irritation at the sensitive nasal bridge area

### Circular diffused venting

directs air away for less user and partner disturbance

### Quick release elbow

with squeeze-tabs to make attaching and detaching quicker and easier

### Flex-wing forehead support


allows a dynamic range of fit and adjustment, and feels light on the face

### Dual wall Spring Air™ cushion

delivers a high-performance seal with even more comfort

### Headgear loops

simplify attachment for improved usability



## Product Codes

Quattro Air mask system Small	62701
Quattro Air mask system Medium	62702
Quattro Air mask system Large	62703
Quattro Air for Her mask system X-Small	62740
Quattro Air for Her mask system Small	62741
Quattro Air for Her mask system Medium	62742

**Quattro™ Air**  
FULL FACE MASK

**Quattro™ Air for Her**  
FULL FACE MASK



**ResMed Ltd**  
1 Elizabeth Macarthur Drive  
Bella Vista NSW 2153 Australia

**Distributed by:** ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. **EC/REP** ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See [www.resmed.com](http://www.resmed.com) for other ResMed locations worldwide. Quattro, SoftEdge and Spring Air are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. © 2013 ResMed Ltd. Specifications may change without notice. 1015994/1 2013-03